

Antibodies Product Datasheet

Product Name:	p73 Antibody
Catalog Number:	D12177N
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human p73 around the non-acetylation site of K327.
Concentration:	1 mg/ml
Gene Name:	TP73
Protein Name:	Tumor protein p73
Other Name:	TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein
Human Swiss Prot No.:	O15350
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	p73 Polyclonal Antibody detects endogenous levels of p73 protein.
Storage and Stability:	-20°C/1 year