

Antibodies Product Datasheet

Product Name:	Dok-2 (Phospho-Tyr345) Antibody
Catalog Number:	A8467
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human Dok-2 (phospho Tyr345)
Concentration:	1 mg/ml
Gene Name:	DOK2
Protein Name:	Docking protein 2
Other Name:	DOK2; Docking protein 2; Downstream of tyrosine kinase 2; p56(dok-2)
Human Swiss Prot No.:	O60496
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Dok-2 (Y345) Polyclonal Antibody detects endogenous levels of Dok-2 protein only when phosphorylated at Y345.
Storage and Stability:	-20°C/1 year