

Antibodies Product Datasheet

Product Name:	PROC (light chain, Cleaved-Leu179) Antibody
Catalog Number:	L0455
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human PROC (light chain, Cleaved-Leu179)
Concentration:	1 mg/ml
Gene Name:	PROC
Protein Name:	PROC (light chain, Cleaved-Leu179)
Other Name:	Vitamin K-dependent protein C (EC 3.4.21.69) (Anticoagulant protein C) (Autoprothrombin IIA) (Blood coagulation factor XIV) [Cleaved into: Vitamin K-dependent protein C light chain; Vitamin K-dependent protein C heavy chain; Activation peptide]
Human Swiss Prot No.:	P04070
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human PROC (light chain, Cleaved-Leu179, protein was cleaved amino acid sequence between 179-180)
Storage and Stability:	-20°C/1 year