

Antibodies Product Datasheet

Product Name:	INT7 Antibody
Catalog Number:	C44267
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human INT7
Concentration:	1 mg/ml
Gene Name:	INTS7 C1orf73
Protein Name:	INT7
Other Name:	Integrator complex subunit 7 (Int7)
Human Swiss Prot No.:	Q9NVH2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.258% sodium azide.
Specificity:	This antibody detects endogenous levels of INT7 at Human/Mouse
Storage and Stability:	-20°C/1 year