

Antibodies Product Datasheet

Product Name:	NFAT3 Antibody
Catalog Number:	R12-3090
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	NFAT3 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NFAT3.
Concentration:	1 mg/ml
Gene Name:	NFATC4
Other Name:	NF-AT3; NF-ATc4; NFAC4; NFATC4; NFC4; Nuclear factor of activated T-cells, cytoplasmic 4; T cell transcription factor NFAT3; Transcription complex subunit NF-ATc4
Human Swiss Prot No.:	Q14934
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	NFAT3 Antibody antibody detects endogenous levels of NFAT3.
Storage and Stability:	Stable at -20°C for at least 1 year.