
Antibodies Product Datasheet

Product Name:	GPR152 Antibody
Catalog Number:	G303
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR152. AA range:301-350
Concentration:	1 mg/ml
Gene Name:	GPR152
Protein Name:	Probable G-protein coupled receptor 152
Other Name:	GPR152; PGR5; Probable G-protein coupled receptor 152; G-protein coupled receptor PGR5
Human Swiss Prot No.:	Q8TDT2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR152 Polyclonal Antibody detects endogenous levels of GPR152 protein.
Storage and Stability:	-20°C/1 year