

Antibodies Product Datasheet

Product Name:	Myosin VI Antibody
Catalog Number:	C16791
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Myosin VI . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	MYO6
Protein Name:	Unconventional myosin-VI
Other Name:	MYO6; KIAA0389; Unconventional myosin-VI; Unconventional myosin-6
Human Swiss Prot No.:	Q9UM54
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Myosin VI Polyclonal Antibody detects endogenous levels of Myosin VI protein.
Storage and Stability:	-20°C/1 year