

Antibodies Product Datasheet

Product Name:	CARP Antibody
Catalog Number:	C11072
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ANKRD1. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	ANKRD1
Protein Name:	Ankyrin repeat domain-containing protein 1
Other Name:	ANKRD1; C193; CARP; HA1A2; Ankyrin repeat domain-containing protein 1; Cardiac ankyrin repeat protein; Cytokine-inducible gene C-193 protein; Cytokine-inducible nuclear protein
Human Swiss Prot No.:	Q15327
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CARP Polyclonal Antibody detects endogenous levels of CARP protein.
Storage and Stability:	-20°C/1 year