

Antibodies Product Datasheet

Product Name:	Angptl7 Antibody
Catalog Number:	C14484
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ANGPTL7. AA range:297-346
Concentration:	1 mg/ml
Gene Name:	ANGPTL7
Protein Name:	Angiopoietin-related protein 7
Other Name:	ANGPTL7; CDT6; Angiopoietin-related protein 7; Angiopoietin-like factor; Angiopoietin-like protein 7; Cornea-derived transcript 6 protein
Human Swiss Prot No.:	O43827
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Angptl7 Polyclonal Antibody detects endogenous levels of Angptl7 protein.
Storage and Stability:	-20°C/1 year