

Antibodies Product Datasheet

Product Name:	MMP24 Antibody
Catalog Number:	C44779
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 520-600
Concentration:	1 mg/ml
Gene Name:	MMP24 MT5MMP
Protein Name:	Matrix metalloproteinase-24 (MMP-24) (EC 3.4.24.-) (Membrane-type matrix metalloproteinase 5) (MT-MMP 5) (MTMMP5) (Membrane-type-5 matrix metalloproteinase) (MT5-MMP) (MT5MMP) [Cleaved into: Processed matrix metalloproteinase-24]
Other Name:	
Human Swiss Prot No.:	Q9Y5R2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MMP24 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year