

Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: TWEAK Antibody

Catalog Number: C10105

Reactivity: Human, Mouse

Applications: WB, IHC-p, IF/ICC, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in

other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human TNFSF12. AA range:41-90

Concentration: 1 mg/ml

Gene Name: TNFSF12

Protein Name: Tumor necrosis factor ligand superfamily member 12

Other Name: TNFSF12; APO3L; DR3LG; Tumor necrosis factor ligand superfamily

member 12; APO3 ligand; TNF-related weak inducer of apoptosis;

TWEAK

Human Swiss Prot No.: O43508

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: TWEAK Polyclonal Antibody detects endogenous levels of TWEAK

protein.

Storage and Stability: -20°C/1 year