

Antibodies Product Datasheet

Product Name:	NANS Antibody
Catalog Number:	C50867
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human NANS
Concentration:	1 mg/ml
Gene Name:	NANS SAS
Protein Name:	Sialic acid synthase (N-acetylneuraminate synthase) (EC 2.5.1.56) (N-acetylneuraminate-9-phosphate synthase) (EC 2.5.1.57) (N-acetylneuraminic acid phosphate synthase) (N-acetylneuraminic acid synthase)
Other Name:	
Human Swiss Prot No.:	Q9NR45
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of NANS at Human
Storage and Stability:	Store at -20° C/1 year