

## Antibodies Product Datasheet

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | PLEK2 Antibody  |
| <b>Catalog Number:</b>        | C45569  |
| <b>Reactivity:</b>            | Human, Mouse  |
| <b>Applications:</b>          | WB  |
| <b>Dilutions:</b>             | WB 1:500-2000   |
| <b>Modification:</b>          | Unmodified/Total  |
| <b>Source:</b>                | Polyclonal Rabbit   |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from human PLEK2  |
| <b>Concentration:</b>         | 1 mg/ml   |
| <b>Gene Name:</b>             | PLEK2   |
| <b>Protein Name:</b>          | PLEK2   |
| <b>Other Name:</b>            | Pleckstrin-2  |
| <b>Human Swiss Prot No.:</b>  | Q9NYT0  |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.38% sodium azide.                                   |
| <b>Specificity:</b>           | This antibody detects endogenous levels of PLEK2 at Human/Mouse   |
| <b>Storage and Stability:</b> | -20°C/1 year  |