

Antibodies Product Datasheet

Product Name:	JAK2 (Phospho-Tyr931) Antibody
Catalog Number:	P12-1129
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	JAK2 (Phospho-Tyr931) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human JAK2 around the phosphorylation site of Tyr931.
Concentration:	1 mg/ml
Gene Name:	JAK2
Other Name:	JAK-2, JAK2, Janus kinase 2, kinase Jak2
Human Swiss Prot No.:	O60674
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	JAK2 (Phospho-Tyr931) Antibody detects endogenous levels of JAK2 only when phosphorylated at Tyr931.
Storage and Stability:	Stable at -20°C for at least 1 year.