
Antibodies Product Datasheet

Product Name:	Anti-Progesterone Receptor (PR) [ABT-PR]
Catalog Number:	V0114
Clone ID:	ABT-PR
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	P06401
Gene Name:	PGR
Other Name:	NR3C3; nuclear receptor 3C3; PGR; PRGR; PR
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-Progesterone Receptor (PR) [ABT-PR] detects endogenous levels of total Progesterone Receptor (PR) protein.
Storage and Stability:	-20°C/1 year