

Antibodies Product Datasheet

Product Name:	IGFBP6 Antibody
Catalog Number:	C30055
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IGFBP6. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	IGFBP6
Protein Name:	Insulin-like growth factor-binding protein 6
Other Name:	IGFBP6; IBP6; Insulin-like growth factor-binding protein 6; IBP-6; IGF-binding protein 6; IGFBP-6
Human Swiss Prot No.:	P24592
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IGFBP6 Polyclonal Antibody detects endogenous levels of IGFBP6 protein.
Storage and Stability:	-20°C/1 year