

Antibodies Product Datasheet

Product Name:	NRAM2 Antibody
Catalog Number:	C45070
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NRAM2
Concentration:	1 mg/ml
Gene Name:	SLC11A2 DCT1 DMT1 NRAMP2 OK/SW-cl.20
Protein Name:	NRAM2
Other Name:	Natural resistance-associated macrophage protein 2 (NRAMP 2) (Divalent cation transporter 1) (Divalent metal transporter 1) (DMT-1)
Human Swiss Prot No.:	P49281
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.156% sodium azide.
Specificity:	This antibody detects endogenous levels of NRAM2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year