

Antibodies Product Datasheet

Product Name:	USP38 Antibody
Catalog Number:	C19383
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human USP38. AA range:551-600
Concentration:	1 mg/ml
Gene Name:	USP38
Protein Name:	Ubiquitin carboxyl-terminal hydrolase 38
Other Name:	USP38; KIAA1891; Ubiquitin carboxyl-terminal hydrolase 38; Deubiquitinating enzyme 38; HP43.8KD; Ubiquitin thioesterase 38; Ubiquitin-specific-processing protease 38
Human Swiss Prot No.:	Q8NB14
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	USP38 Polyclonal Antibody detects endogenous levels of USP38 protein.
Storage and Stability:	-20°C/1 year