

Antibodies Product Datasheet

Product Name:	Cystatin C Monoclonal Antibody
Catalog Number:	N0178
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Cystatin C expressed in E. Coli.
Concentration:	
Gene Name:	CST3
Protein Name:	Cystatin-C
Other Name:	CST3; Cystatin-C; Cystatin-3; Gamma-trace; Neuroendocrine basic polypeptide; Post-gamma-globulin
Human Swiss Prot No.:	P01034
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Cystatin C Monoclonal Antibody detects endogenous levels of Cystatin C protein.
Storage and Stability:	-20°C/1 year