

Antibodies Product Datasheet

Product Name:	CNOT7 Antibody
Catalog Number:	C11858
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CNOT7. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	CNOT7
Protein Name:	CCR4-NOT transcription complex subunit 7
Other Name:	CNOT7; CAF1; CCR4-NOT transcription complex subunit 7; BTG1-binding factor 1; CCR4-associated factor 1; CAF-1
Human Swiss Prot No.:	Q9UIV1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CNOT7 Polyclonal Antibody detects endogenous levels of CNOT7 protein.
Storage and Stability:	-20°C/1 year