

Antibodies Product Datasheet

Product Name:	CDX2 mouse mAb
Catalog Number:	N0912
Reactivity:	Human,Monkey
Applications:	WB,IP,ICC
Dilutions:	wb 1:1000 icc 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human CDX2 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	cdx2
Protein Name:	
Other Name:	Caudal type homeo box 2;Caudal type homeo box transcription factor 2;Caudal type homeobox 2;Caudal type homeobox protein 2;caudal type homeobox transcription factor 2; Caudal-type homeobox protein 2;CDX 2;CDX 3;CDX-3;Cdx2;CDX2_HUMAN;CDX3;Homeobox protein
Human Swiss Prot No.:	Q99626
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CDX2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year