

Antibodies Product Datasheet

Product Name:	KRT37/38 Antibody
Catalog Number:	R12-2958
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	KRT37/38 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KRT37/38.
Concentration:	1 mg/ml
Gene Name:	KRT37/KRT38
Other Name:	Keratin type I cuticular Ha7; Hair keratin type I Ha7; Keratin-37; K37; KRT37; HHA7; HKA7; KRTHA7; Keratin type I cuticular Ha8; Hair keratin type I Ha8; Keratin-38; KRT38; HHA8; HKA8; KRTHA8
Human Swiss Prot No.:	O76014/O76015
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	KRT37/38 Antibody antibody detects endogenous levels of KRT37/38.
Storage and Stability:	Stable at -20°C for at least 1 year.