

Antibodies Product Datasheet

Product Name:	ARHGAP12 Antibody
Catalog Number:	C18366
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ARHGAP12 . at AA range: 310-390
Concentration:	1 mg/ml
Gene Name:	ARHGAP12
Protein Name:	Rho GTPase-activating protein 12
Other Name:	ARHGAP12; Rho GTPase-activating protein 12; Rho-type GTPase-activating protein 12
Human Swiss Prot No.:	Q8IWW6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARHGAP12 Polyclonal Antibody detects endogenous levels of ARHGAP12 protein.
Storage and Stability:	-20°C/1 year