

Antibodies Product Datasheet

Product Name:	Krs-1/2 Antibody
Catalog Number:	B0687
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Mst1/2. AA range:149-198
Concentration:	1 mg/ml
Gene Name:	STK3/STK4
Protein Name:	Serine/threonine-protein kinase 3/4
Other Name:	STK3; KRS1; MST2; Serine/threonine-protein kinase 3; Mammalian STE20-like protein kinase 2; MST-2; STE20-like kinase MST2; Serine/threonine-protein kinase Krs-1; STK4; KRS2; MST1; Serine/threonine-protein kinase 4; Mammalian STE20-like prot
Human Swiss Prot No.:	Q13188/Q13043
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Krs-1/2 Polyclonal Antibody detects endogenous levels of Krs-1/2 protein.
Storage and Stability:	-20°C/1 year