
Antibodies Product Datasheet

Product Name:	PHKB Antibody
Catalog Number:	C11458
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KPBB. AA range:661-710
Concentration:	1 mg/ml
Gene Name:	PHKB
Protein Name:	Phosphorylase b kinase regulatory subunit beta
Other Name:	PHKB; Phosphorylase b kinase regulatory subunit beta; Phosphorylase kinase subunit beta
Human Swiss Prot No.:	Q93100
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PHKB Polyclonal Antibody detects endogenous levels of PHKB protein.
Storage and Stability:	-20°C/1 year