
Antibodies Product Datasheet

Product Name:	Tubulin β mouse mAb(PT0354)
Catalog Number:	N1728
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),WB
Dilutions:	IHC-p 1:100-300, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Tubulin β
Concentration:	1 mg/ml
Gene Name:	TUBB3 TUBB4
Protein Name:	Tubulin β
Other Name:	Tubulin beta-3 chain (Tubulin beta-4 chain) (Tubulin beta-III)
Human Swiss Prot No.:	Q13509
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Tubulin β
Storage and Stability:	-20°C/1 year