

Antibodies Product Datasheet

Product Name:	SLC38A1 Antibody
Catalog Number:	C48297
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SLC38A1
Concentration:	1 mg/ml
Gene Name:	SLC38A1 ATA1 NAT2 SAT1 SNAT1
Protein Name:	Sodium-coupled neutral amino acid transporter 1 (Amino acid transporter A1) (N-system amino acid transporter 2) (Solute carrier family 38 member 1) (System A amino acid transporter 1) (System N amino acid transporter 1)
Other Name:	
Human Swiss Prot No.:	Q9H2H9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SLC38A1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year