

Antibodies Product Datasheet

Product Name:	GAL3ST4 Antibody
Catalog Number:	C50602
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GAL3ST4
Concentration:	1 mg/ml
Gene Name:	GAL3ST4 PP6968
Protein Name:	Galactose-3-O-sulfotransferase 4 (Gal3ST-4) (EC 2.8.2.-) (Beta-galactose-3-O-sulfotransferase 4) (Gal-beta-1,3-GalNAc 3'-sulfotransferase)
Other Name:	
Human Swiss Prot No.:	Q96RP7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GAL3ST4 at Human
Storage and Stability:	Store at -20° C/1 year