

## Antibodies Product Datasheet

<b>Product Name:</b>	PTEN (Phospho-Ser380) Antibody
<b>Catalog Number:</b>	A7202
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human PTEN around the phosphorylation site of Ser380. AA range:366-395
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PTEN
<b>Protein Name:</b>	Phosphatidylinositol 3,4,5-trisphosphate 3-phosphatase and dual-specificity protein phosphatase PTEN
<b>Other Name:</b>	PTEN; MMAC1; TEP1; Phosphatidylinositol 3; 4,5-trisphosphate 3-phosphatase and dual-specificity protein phosphatase PTEN; Mutated in multiple advanced cancers 1; Phosphatase and tensin homolog
<b>Human Swiss Prot No.:</b>	P60484
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-PTEN (S380) Polyclonal Antibody detects endogenous levels of PTEN protein only when phosphorylated at S380.
<b>Storage and Stability:</b>	-20°C/1 year