

Antibodies Product Datasheet

Product Name:	CD81 Antibody
Catalog Number:	C30385
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CD81. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	CD81
Protein Name:	CD81 antigen
Other Name:	CD81; TAPA1; TSPAN28; CD81 antigen; 26 kDa cell surface protein TAPA-1; Target of the antiproliferative antibody 1; Tetraspanin-28; Tspan-28; CD81
Human Swiss Prot No.:	P60033
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD81 Polyclonal Antibody detects endogenous levels of CD81 protein.
Storage and Stability:	-20°C/1 year