

Antibodies Product Datasheet

Product Name:	RPGF4 Antibody
Catalog Number:	C46000
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 240-320
Concentration:	1 mg/ml
Gene Name:	RAPGEF4 CGEF2 EPAC2
Protein Name:	Rap guanine nucleotide exchange factor 4 (Exchange factor directly activated by cAMP 2) (Exchange protein directly activated by cAMP 2) (EPAC 2) (cAMP-regulated guanine nucleotide exchange factor II) (cAMP-GEFII)
Other Name:	
Human Swiss Prot No.:	Q8WZA2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	RPGF4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year