
Antibodies Product Datasheet

Product Name:	Cleaved-COL3A1 (G1221) Antibody
Catalog Number:	L0211
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Collagen III alpha1. AA range:1172-1221
Concentration:	1 mg/ml
Gene Name:	COL3A1
Protein Name:	Collagen alpha-1(III) chain
Other Name:	COL3A1; Collagen alpha-1(III) chain
Human Swiss Prot No.:	P02461
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-COL3A1 (G1221) Polyclonal Antibody detects endogenous levels of fragment of activated COL3A1 protein resulting from cleavage adjacent to G1221.
Storage and Stability:	-20°C/1 year