

Antibodies Product Datasheet

Product Name:	DLG2 Antibody
Catalog Number:	R12-2701
Reactivity:	Human
Applications:	WB, IHC, ELISA
Dilutions:	DLG2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DLG2.
Concentration:	1 mg/ml
Gene Name:	DLG2
Other Name:	channel associated protein of synapse-110; chapsyn-110; Discs large 2; DLG2; postsynaptic density protein PSD-93
Human Swiss Prot No.:	Q15700
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	DLG2 Antibody antibody detects endogenous levels of DLG2.
Storage and Stability:	Stable at -20°C for at least 1 year.