

Antibodies Product Datasheet

Product Name:	RPB11 Antibody
Catalog Number:	R12-3448
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	RPB11 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RPB11.
Concentration:	1 mg/ml
Gene Name:	POLR2J
Other Name:	DNA-directed RNA polymerase II subunit RPB11-a; RNA polymerase II subunit B11-a; RPB11a; DNA-directed RNA polymerase II subunit J-1; RNA polymerase II 13.3 kDa subunit; POLR2J; POLR2J1
Human Swiss Prot No.:	P52435
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	RPB11 Antibody antibody detects endogenous levels of RPB11.
Storage and Stability:	Stable at -20°C for at least 1 year.