

Antibodies Product Datasheet

Product Name:	c-Fms Antibody
Catalog Number:	B0876
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CSFR. AA range:781-830
Concentration:	1 mg/ml
Gene Name:	CSF1R
Protein Name:	Macrophage colony-stimulating factor 1 receptor
Other Name:	CSF1R; FMS; Macrophage colony-stimulating factor 1 receptor; CSF-1 receptor; CSF-1-R; CSF-1R; M-CSF-R; Proto-oncogene c-Fms; CD antigen CD115
Human Swiss Prot No.:	P07333
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	c-Fms Polyclonal Antibody detects endogenous levels of c-Fms protein.
Storage and Stability:	-20°C/1 year