

Antibodies Product Datasheet

Product Name:	Smad4 Antibody
Catalog Number:	R12-3546
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Smad4 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Smad4.
Concentration:	1 mg/ml
Gene Name:	SMAD4
Other Name:	Mothers against decapentaplegic homolog 4; MAD homolog 4; Mothers against DPP homolog 4; Deletion target in pancreatic carcinoma 4; SMAD family member 4; SMAD 4; Smad4; hSMAD4; SMAD4; DPC4; MADH4
Human Swiss Prot No.:	Q13485
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Smad4 Antibody antibody detects endogenous levels of Smad4.
Storage and Stability:	Stable at -20°C for at least 1 year.