
Antibodies Product Datasheet

Product Name:	Cdc5L Antibody
Catalog Number:	C12171
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDC5L. AA range:721-770
Concentration:	1 mg/ml
Gene Name:	CDC5L
Protein Name:	Cell division cycle 5-like protein
Other Name:	CDC5L; KIAA0432; PCDC5RP; Cell division cycle 5-like protein; Cdc5-like protein; Pombe cdc5-related protein
Human Swiss Prot No.:	Q99459
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc5L Polyclonal Antibody detects endogenous levels of Cdc5L protein.
Storage and Stability:	-20°C/1 year