

## Antibodies Product Datasheet

<b>Product Name:</b>	IL-9R (Phospho-Ser519) Antibody
<b>Catalog Number:</b>	A1067
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, IHC-p, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human IL-9R around the phosphorylation site of Ser519. AA range:472-521
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	IL9R
<b>Protein Name:</b>	Interleukin-9 receptor
<b>Other Name:</b>	IL9R; Interleukin-9 receptor; IL-9 receptor; IL-9R; CD antigen CD129
<b>Human Swiss Prot No.:</b>	Q01113
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-IL-9R (S519) Polyclonal Antibody detects endogenous levels of IL-9R protein only when phosphorylated at S519.
<b>Storage and Stability:</b>	-20°C/1 year