

Antibodies Product Datasheet

Product Name:	CISD2 Antibody
Catalog Number:	C43097
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 31-80
Concentration:	1 mg/ml
Gene Name:	CISD2 CDGSH2 ERIS ZCD2
Protein Name:	CDGSH iron-sulfur domain-containing protein 2 (Endoplasmic reticulum intermembrane small protein) (MitoNEET-related 1 protein) (Miner1) (Nutrient-deprivation autophagy factor-1) (NAF-1)
Other Name:	
Human Swiss Prot No.:	Q8N5K1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CISD2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year