

Antibodies Product Datasheet

Product Name:	AKAP 79 Antibody
Catalog Number:	C10472
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AKAP5. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	AKAP5
Protein Name:	A-kinase anchor protein 5
Other Name:	AKAP5; AKAP79; A-kinase anchor protein 5; AKAP-5; A-kinase anchor protein 79 kDa; AKAP 79; H21; cAMP-dependent protein kinase regulatory subunit II high affinity-binding protein
Human Swiss Prot No.:	P24588
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AKAP 79 Polyclonal Antibody detects endogenous levels of AKAP 79 protein.
Storage and Stability:	-20°C/1 year