

Antibodies Product Datasheet

Product Name:	Cytochrome P450 2D6 Antibody
Catalog Number:	R12-2678
Reactivity:	Human
Applications:	WB, IHC, ELISA
Dilutions:	Cytochrome P450 2D6 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 2D6.
Concentration:	1 mg/ml
Gene Name:	CYP2D6
Other Name:	Cytochrome P450 2D6; CYP11D6; Cytochrome P450-DB1; Debrisoquine 4-hydroxylase; CYP2D6; CYP2DL1
Human Swiss Prot No.:	P10635
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Cytochrome P450 2D6 Antibody antibody detects endogenous levels of Cytochrome P450 2D6.
Storage and Stability:	Stable at -20°C for at least 1 year.