

Antibodies Product Datasheet

Product Name:	MDM2 Antibody
Catalog Number:	B0508
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MDM2. AA range:132-181
Concentration:	1 mg/ml
Gene Name:	MDM2
Protein Name:	E3 ubiquitin-protein ligase Mdm2
Other Name:	MDM2; E3 ubiquitin-protein ligase Mdm2; Double minute 2 protein; Hdm2; Oncoprotein Mdm2; p53-binding protein Mdm2
Human Swiss Prot No.:	Q00987
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MDM2 Polyclonal Antibody detects endogenous levels of MDM2 protein.
Storage and Stability:	-20°C/1 year