

Antibodies Product Datasheet

Product Name:	BM88 Antibody
Catalog Number:	C12185
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CEND. AA range:11-60
Concentration:	1 mg/ml
Gene Name:	CEND1
Protein Name:	Cell cycle exit and neuronal differentiation protein 1
Other Name:	CEND1; BM88; Cell cycle exit and neuronal differentiation protein 1; BM88 antigen
Human Swiss Prot No.:	Q8N111
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BM88 Polyclonal Antibody detects endogenous levels of BM88 protein.
Storage and Stability:	-20°C/1 year