

Antibodies Product Datasheet

Product Name:	TIP60 Antibody
Catalog Number:	C11451
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TIP60. AA range:371-420
Concentration:	1 mg/ml
Gene Name:	KAT5
Protein Name:	Histone acetyltransferase KAT5
Other Name:	KAT5; HTATIP; TIP60; Histone acetyltransferase KAT5; 60 kDa Tat-interactive protein; Tip60; Histone acetyltransferase HTATIP; HIV-1 Tat interactive protein; Lysine acetyltransferase 5; cPLA(2)-interacting protein
Human Swiss Prot No.:	Q92993
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TIP60 Polyclonal Antibody detects endogenous levels of TIP60 protein.
Storage and Stability:	-20°C/1 year