
Antibodies Product Datasheet

Product Name:	TIMP3 Antibody
Catalog Number:	R12-3625
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	TIMP3 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TIMP3.
Concentration:	1 mg/ml
Gene Name:	TIMP3
Other Name:	Metalloproteinase inhibitor 3; Tissue inhibitor of metalloproteinases 3; TIMP-3; Protein MIG-5
Human Swiss Prot No.:	P35625
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	TIMP3 Antibodyantibody detects endogenous levels of TIMP3.
Storage and Stability:	Stable at -20°C for at least 1 year.