

Antibodies Product Datasheet

Product Name:	Anti-Human Mammaglobin (ABT079) Mouse Monoclonal Antibody
Catalog Number:	V0389
Clone ID:	ABT079
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	SCGB2A2
Gene Name:	Q13296
Other Name:	Mammaglobin 1, Mammaglobin-1, Mammaglobin-A, MGB1, Scgb2a2, Secretoglobin family 2A member 2, SG2A2, SG2A2_HUMAN, UGB2
Isotype	IgG1, Kappa
Localization:	Breast hyperplasia
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Mammaglobin protein.
Storage and Stability:	-20°C/1 year