

Antibodies Product Datasheet

Product Name:	OATP1 Antibody
Catalog Number:	C18884
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLCO1A2. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	SLCO1A2
Protein Name:	Solute carrier organic anion transporter family member 1A2
Other Name:	SLCO1A2; OATP; OATP1; OATP1A2; SLC21A3; Solute carrier organic anion transporter family member 1A2; OATP-A; Organic anion-transporting polypeptide 1; OATP-1; Sodium-independent organic anion transporter; Solute carrier family 21 member 3
Human Swiss Prot No.:	P46721
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	OATP1 Polyclonal Antibody detects endogenous levels of OATP1 protein.
Storage and Stability:	-20°C/1 year