

## Antibodies Product Datasheet

|                               |  |
|-------------------------------|--|
| <b>Product Name:</b>          | TASP1 Antibody   |
| <b>Catalog Number:</b>        | C49038   |
| <b>Reactivity:</b>            | Human, Mouse   |
| <b>Applications:</b>          | WB   |
| <b>Dilutions:</b>             | WB 1:500-2000  |
| <b>Modification:</b>          | Unmodified/Total   |
| <b>Source:</b>                | Polyclonal, Rabbit,IgG   |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                |
| <b>Immunogen:</b>             | Synthesized peptide derived from human TASP1   |
| <b>Concentration:</b>         | 1 mg/ml  |
| <b>Gene Name:</b>             | TASP1 C20orf13   |
| <b>Protein Name:</b>          | Threonine aspartase 1 (Taspase-1) (EC 3.4.25.-) [Cleaved into:<br>Threonine aspartase subunit alphaThreonine aspartase subunit beta] |
| <b>Other Name:</b>            |  |
| <b>Human Swiss Prot No.:</b>  | Q9H6P5   |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.  |
| <b>Specificity:</b>           | This antibody detects endogenous levels of TASP1 at Human, Mouse   |
| <b>Storage and Stability:</b> | Store at -20° C/1 year   |