

Antibodies Product Datasheet

Product Name:	PIP5KIII Antibody
Catalog Number:	C11935
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PIP5K. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	PIKFYVE
Protein Name:	1-phosphatidylinositol 3-phosphate 5-kinase
Other Name:	PIKFYVE; KIAA0981; PIP5K3; 1-phosphatidylinositol 3-phosphate 5-kinase; Phosphatidylinositol 3-phosphate 5-kinase; FYVE finger-containing phosphoinositide kinase; PIKfyve; Phosphatidylinositol 3-phosphate 5-kinase type III; PIPkin-III; Type
Human Swiss Prot No.:	Q9Y2I7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PIP5KIII Polyclonal Antibody detects endogenous levels of PIP5KIII protein.
Storage and Stability:	-20°C/1 year