

Antibodies Product Datasheet

Product Name:	MYH1 Antibody
Catalog Number:	C44870
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	MYH1
Protein Name:	Myosin-1 (Myosin heavy chain 1) (Myosin heavy chain 2x) (MyHC-2x) (Myosin heavy chain IIx/d) (MyHC-IIx/d) (Myosin heavy chain, skeletal muscle, adult 1)
Other Name:	
Human Swiss Prot No.:	P12882
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MYH1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year